

## Protocol: Cell Cycle/Whole Cell

### Reagents:

- 1.) Ice cold 100% Ethanol
- 2.) Propidium Iodide (PI) staining solution: To 10ml of PBS, add 500 $\mu$ l of a 1mg/ml solution of PI, 2mg of RNase A, and 100 $\mu$ l of 10% Triton-X solution (final concentrations: 50 $\mu$ g/ml PI, 0.1% Triton-X).

### Procedure:

#### Ethanol fixation:

1. Wash  $10^6$  cells twice in 2ml PBS. (Centrifuge at  $\cong 200 \times g$  for 5 minutes, decant, vortex, repeat).
2. Thoroughly resuspend cells in 200 $\mu$ l PBS-it is essential that cells be in single cell suspension prior to fixation with ethanol.
3. Add 800 $\mu$ l of ice-cold 100% ethanol drop-wise while gentle vortexing.
4. Keep cells in fixative  $\geq 2$  hours.

**Note:** Cells suspended in 80% ethanol can be stored at  $-20^{\circ}\text{C}$  for several weeks.

#### PI staining: (This step can cause lots of cell loss)

5. Add 1ml of PBS to cells.
6. Spin at 800g for 5min
7. Aspirate and repeat wash with 2ml PBS
8. Resuspend cell pellet in 500 $\mu$ l PI staining solution; vortex and incubate for 20minutes at  $37^{\circ}\text{C}$ .
9. Cover with foil and store at  $4^{\circ}\text{C}$  until ready to analyze on flow cytometer.