

Protocol: Cell Cycle/Whole Cell

Reagents:

- 1.) Ice cold 100% Ethanol
- 2.) Propidium Iodide (PI) staining solution: To 10ml of PBS, add 500ul of a 1mg/ml solution of PI, 2mg of RNase A, and 100ul of 10% Triton-X solution (final concentrations: 50ug/ml PI, 0.1% Triton-X).

Procedure:

Ethanol fixation:

1. Wash 10^6 cells twice in 2ml PBS. (Centrifuge at $\approx 300 \times g$ for 5 minutes, decant, vortex, repeat).
2. Thoroughly resuspend cells in 200ul PBS-it essential that cells be in single cell suspension prior to fixation with ethanol.
3. Add 800ul of ice-cold 100% ethanol dropwise while gentle vortexing.
4. Keep cells in fixative ≥ 2 hours.

Note: Cells suspended in 80% ethanol can be stored at -20°C for several weeks.

PI staining: (This step can cause lots of cell loss)

5. Add 1ml of PBS to cells.
6. Spin at 800g for 5min
7. Aspirate, and repeat wash with 2ml PBS
8. Resuspend cell pellet in 500ml PI staining solution; vortex and incubate for 20minutes at 37°C .
9. Cover with foil and store at 4°C until ready to analyze on flow cytometer.