

Cytometer	HCI Aria					
Software	BD FACSDiva 8.0					
Serial #	P69500142					
Laser	Wavelength	Filter	Parameter Names			
Blue	488nm	488/10	SSC			
		525/50	FITC	GFP		B (525/50)
		695/40	PerCP	PerCP Cy5.5		B (695/40)
Yel-Gr	561nm	582/15	PE			YG (582/15)
		610/20	PE Texas Red	PE eFluor 610	PE-CF594	YG (610/20)
		660/20	PE Cy5			YG (660/20)
		710/50	PE Cy5.5			YG (710/50)
		780/60	PE Cy7			YG (780/60)
Red	640nm	660/20	APC			R (660/20)
		730/45	AF700			R (730/45)
		780/60	APC Cy7	APC eFluor 780	APC H7	R (780/60)
Violet	408nm	450/50	BV421			V (450/50)
		525/50	BV510			V (525/50)
		610/20	BV605			V (610/20)
		670/30	BV650			V (670/30)
		710/50	BV711			V (710/50)
		780/40	BV786			V (780/40)
UV	350nm	380/30	BUV395			UV (380/30)
		515/20	BUV496			UV (515/20)

Cytometer	Avalon	
Block-1		
488nm	525/50	GFP
561nm	585/15	PE
	610/20	RFP
	660LP	Percp
Block-2		
488nm	525/50	GFP
561nm	593/40	PE
	660/40	APC
	780/60	PECy-7

Cytometer	Wintrobe Aria					
Software	BD FACSDiva 8.0					
Serial #	P28200193					
Laser	Wavelength	Filter	Parameter Names			
Blue	488nm	488/10	SSC			
		530/30	FITC	GFP		B (530/30)
		695/40	PerCP	PerCP Cy5.5		B (695/40)
Yel-Gr	561nm	590/20	PE			YG (590/20)
		615/25	PE Texas Red	PE eFluor 610	PE-CF594	YG (615/25)
		667/30	PE Cy5			YG (667/30)
		780/60	PE Cy7			YG (780/60)
Red	640nm	660/20	APC			R (660/20)
		730/45	AF700			R (730/45)
		780/60	APC Cy7	APC eFluor 780	APC H7	R (780/60)
Violet	408nm	450/50	BV421	DAPI		V (450/50)
		525/50	BV510			V (525/50)
		610/20	BV605			V (610/20)
		670/30	BV650			V (670/30)
		710/50	BV711			V (710/50)
		780/60	BV786			V (780/40)

Cytometer	HCI Fortessa			
Software	BD FACSDiva 8.0.1			
Serial #	H64717700146			
Laser Name	Wavelength	Filter	Parameter Names	
Violet	405	780/60	BV786	
		710/50	BV 711	
		670/30	BV 650	
		610/20	BV 605	
		525/50	BV 510	
		450/50	BV 421	
Yel-Gr	561	780/60	PE-Cy7	
		710/50	PE-Cy5.5	
		670/30	PE-Cy5	
		610/20	PE-Texas Red	
		586/15	PE	
Blue	488	710/50	PerCP-Cy5-5	
		530/30	FITC	
		488/10	SSC	
Red	640	780/60	APC-Cy7	
		730/45	Alexa Fluor 700	
		670/30	APC	
UV	350	520/30	BUV 496	
		379/34	BUV 395	

Cytometer	Wintrobe Canto					
Software	BD FACSDiva 6.1.3					
Serial #	V96300741					
Laser Name	Wavelength	Filter	Parameter Names			
Violet/Red	405/640	780/60 BP	APC-Cy7	APC-H7	R (780/60)	
These lasers are Colinear		670/30 BP	APC	Alexa 647	R (670/30)	
		525/50 BP	AmCyan	Live/Dead Aqua	V (525/50)	
		450/50 BP	Alexa Fluor 405	DAPI	Pacific Blue	V (450/50)
Blue	488	670 LP	PE-Cy5	PerCP-Cy5-5	PI	B (670LP)
		530/30 BP	FITC			
			GFP			
			B (530/30			
			CFSE			
		488/10 BP	SSC			
YG-561	561	780/60 BP	PE-Cy7	YG (780/60)		
		585/15 BP	PE	RFP	YG (586/15)	

Cytometer	Celesta				
Software	BD FACSDiva 8.0.1.1				
Serial #	H66034400048				
Laser Name	Wavelength	Filter	Parameter Names		
Violet	405	780/60 BP	BV786		
		670/30 BP	BV650		
		610/20 BP	BV605		
		525/50 BP	BV510	V500	
		450/40 BP	BV421	V450	Pacific Blue
Blue	488	780/60 BP	PE-Cy7		
		695/40 BP	7-AAD	PerCP-Cy5-5	PI
		575/25 BP	PE		
		530/30 BP	FITC		
		488/10 BP	SSC		
Red	640	780/60 BP	APC-Cy7	APC-H7	
		730/45 BP	APC-R700	Alexa Fluor 700	
		670/30 BP	APC	Alexa Fluor 647	

Cytometer	Cytoflex HCl				
Serial Number	BB09006				
Laser Name	Wavelength	Filter	Parameter		
Red	633	763/43	APC-Cy7		R 763/43
		712/25	Alexa Fluor 700		R 712/25
		660/10	APC		R 660/10
Blue	488	690/50	PerCP	PerCP-Cy5-5	B 690/50
		610/20			B 610/20
		525/40	FITC	GFP	B 525/40
		488/10	SSC		B 488/10
Violet	405	763/43	BV711	BV786	V 763/43
		660/10	BV650		V 660/10
		610/20	BV605		V 610/20
		525/40	BV510		V 525/40
		450/45	BV421	DAPI	V 450/45
Yel-Grn	561	763/43	PE-Cy7		YG 763/43
		710/50	PE-Cy5.5		YG 710/50
		675/30	PE-Cy5		YG 675/30
		610/20	PE-CF594	PE-Texas Red	YG 610/20
		585/42	PE		YG 585/42

Cytometer	DXP					
Serial Number						
Laser Name	Wavelength	Filter	Parameter Names			
Violet	405	525/50 BP	BV510	Live/Dead Aqu	V (525/50)	
		450/50 BP	BV421	DAPI	Pacific Blue	V (450/50)
Blue	488	780/60 BP	PE-Cy7	YG (780/60)		
		670 LP	PerCP-Cy5-5	PI		
		585/15 BP	PE			
		530/30 BP	FITC	GFP		
		488/10 BP	SSC			
Red	640	780/60 BP	APC-Cy7	APC-H7	R (780/60)	
		670/30 BP	APC	Alexa 647	R (670/30)	

Cytometer	Cytoflex Genetics				
Serial Number	ba36265				
Laser Name	Wavelength	Filter	Parameter		
Violet	405	525/50 BP	BV510	Live/Dead Aqua	
		450/50 BP	BV421	DAPI	Pacific Blue
Blue	488	780/60 BP	PE-Cy7		
		670 LP	PerCP-Cy5-5	PI	
		585/15 BP	PE		
		530/30 BP	FITC	GFP	
		488/10 BP	SSC		
Red	640	780/60 BP	APC-Cy7	APC-H7	
		670/30 BP	APC	Alexa 647	

Cytometer	Cytoflex Crocker				
Serial Number	BB21071				
Laser Name	Wavelength	Filter	Parameter Names		
Violet	405	525/40 BP	BV510	Live/Dead Aqua	
		450/45 BP	BV421	DAPI	Pacific Blue
Blue	488	690/50 BP	PerCP-Cy5-5	PI	
		525/40 BP	FITC	GFP	
		488/8 BP	SSC		
YG-561	561	780/60 BP	PE-Cy7	YG (780/60)	
		610/20 BP	PE-TxRed		
		585/42 BP	PE	RFP	
Red	640	780/60 BP	APC-Cy7	APC-H7	
		660/10 BP	APC	Alexa 647	

Cytometer	Cytoflex Crocker				
Serial Number	BB21071				
Laser Name	Wavelength	Filter	Parameter Names		
Violet	405	525/40 BP	BV510	Live/Dead Aqua	
		450/45 BP	BV421	DAPI	Pacific Blue
Blue	488	690/50 BP	PerCP-Cy5-5	PI	
		525/40 BP	FITC	GFP	
		488/8 BP	SSC		
YG-561	561	780/60 BP	PE-Cy7	YG (780/60)	
		610/20 BP	PE-TxRed		
		585/42 BP	PE	RFP	
Red	640	780/60 BP	APC-Cy7	APC-H7	
		660/10 BP	APC	Alexa 647	